

Patient information

How to Fit your Earmould and Position your Ear Level Hearing Aid

Department of Audiology ENT Care Group

It is important to get the earmould in your ear properly. Your audiologist will show you how to do this but it will take practice at home before you are an expert.

You may find it helpful to use a mirror when practising, but not everyone finds this to be the case.



Hold the earmould at the back with your finger and thumb. Pull it back past your ear.



Position the ear mould tip into your ear canal.



Push the ear mould down the ear canal, as well as into the bowl of your ear.



Then put in place the bit than goes into the crease at the top of the ear

Use the other hand to pull down your ear lobe. Push your earmould firmly in.



Correctly fitted aid.



Top part of the earmould not properly tucked into the crease at the top of the ear – incorrect.

Feedback

Your feedback is important to us and helps us influence care in the future.

Following your discharge from hospital or attendance at your outpatient appointment you will receive a text asking if you would recommend our service to others. Please take the time to text back, you will not be charged for the text and can opt out at any point. Your co-operation is greatly appreciated.

Further information

**If you require any further information please contact the Audiology Department
Alexandra Wing
Broadgreen Hospital
Thomas Drive L14 3lb**

Tel: 0151 706 2533

Text phone number: 18001 0151 706 2533

Audiology Text Number: 0782 687 3856

Audiologyoutpatientqueries@LiverpoolFT.Nhs.Uk

All Trust approved information is available on request in alternative formats, including other languages, easy read, large print, audio, Braille, moon and electronically.

يمكن توفير جميع المعلومات المتعلقة بالمرضى الموافق عليهم من قبل انتمان المستشفى عند الطلب بصيغ أخرى، بما في ذلك لغات أخرى وبطرق تسهل قراءتها وبالحروف الطباعية الكبيرة وبالصوت وبطريقة برايل للمكفوفين وبطريقة مون والكترونياً.

所有經信托基金批准的患者資訊均可以其它格式提供，包括其它語言、易讀閱讀軟件、大字

體、音頻、盲文、穆恩體 (Moon) 盲文和電子格式，敬請索取。

در صورت تمایل میتوانید کلیه اطلاعات تصویب شده توسط اتحادیه در رابطه با بیماران را به اشکال مختلف در دسترس داشته باشید، از جمله به زبانهای دیگر، به زبان ساده، چاپ درشت، صوت، خط مخصوص کوران، مون و بصورت روی خطی موجود است.

زانیاری پیوهندیدار بهو نهخوشانهی له لایمن تراستهوه پسهند کراون، نهگس داوا بکریت له فورماتهکانی تردا بریتی له زمانهکانی تر، نیزی رید (هاسان خویندنهوه)، چاپی گهوره، شریتی دهنگ، هیلی موون و نهلیکترونیکی ههیه.

所有经信托基金批准的患者信息均可以其它格式提供，包括其它语言、易读阅读软件、大字体、音频、盲文、穆恩体 (Moon) 盲文和电子格式，敬请索取。

Dhammaan warbixinta bukaanleyda ee Ururka ee la oggol yahay waxaa marka la codsado lagu heli karaa nuskhado kale, sida luqado kale, akhris fudud, far waaweyn, dhegeysi, farta braille ee dadka indhaha la', Moon iyo nidaam eletaroonig ah.